
Apresys®

USER MANUAL

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1. Introduction

Thanks for using our hand-held breath alcohol analyzer: AP2020c. To have a better experience on this product, please go through this user manual in great detail.

AP2020c is a fuel cell sensor based alcohol test with the following features:

- (1) UK made Electrochemical/Fuel cell sensor with high reliability;
- (2) Quick and Precise Analysis;
- (3) Simple structured, of small volume and easy carryover.

Prepare for test
Blowing finished

Blowing
Test result

2. Instrument and Its Accessories / 3. Safety and Precautions

2. Instrument and Its Accessories:

- | | |
|----------------|-------|
| 1、 AP2020c | 1unit |
| 2、 Mouthpieces | 10PCS |


Please contact the supplier if you find anything missed or damaged.

3. Safety and Precautions

Pay attention to all the CAUTION after mark to avoid any misoperation deviation and damage.

Please don't try to open the instrument and do repair without our instruction.


4. Products Pictures and Description

5. Operating Instruction

5.1 Startup and shut down

Startup: hold the ON/OFF button till a “beep” has been given. The device will come to test interface:


Shut down: hold the ON/OFF button till a “beep” has been given.


5.2 Test operation

Remove a disposable mouthpiece from the plastic wrapping, insert into air inlet of AP2020c. Take a deep breath and blow, moderately and continuously into the mouthpiece. The pressure scroll bar is reaching towards the end and beep tone is heard indicating a proper breath flow. The test result will be given:


5. Operating Instruction


Prepare for test

Blowing finished

Blowing

Test result

If the driver stopped blowing before the pressure scroll bar reaching the end, a warning tone is heard and the screen of the device will display “blowing interrupted”, a re-test will be required. If no alcohol detected, LED will be green, If alcohol detected, LED will be Blue or Red.




5. Operating Instruction

CAUTION:

Please make sure the time interval between two continuous tests >2 minute. If the test result is more than BAC 80mg/100mL (BrAC 0.36mg/L) , the next test can not be conducted within 3 minutes.

5.3 Battery low power warning

Please change the battery immediately when low power warning flashing on the screen, the instrument will turn off automatically when the second low power warning occur.

Please change the battery immediately when low power warning flashing on the screen.

6. Technical Data

Details	Technical Data	
Measurement Principle	Electrochemical sensor	
Measurement Range	0.000mg/L ~ 1.500 mg/L	
Resolution	0.001 mg/L	
Accuracy: Measurement Error, in relation to the ethanol standard	Measurement Range C / (mg/L)	Measurement Error
	$C < 0.400$	± 0.020 mg/L
	$0.400 \leq C < 1.000$	$\pm 5\%$
	$C \geq 1.000$	$\pm 20\%$
Calibration	Every 6 months	
Operation temperature	0℃ ~ 40℃	
Storage temperature	-20℃ ~ 70℃	
Optimum Storage Conditions	Temperature: $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ Humidity: <80%RH Pressure: 86 kpa ~ 106 kpa No vapor or electromagnetic disturbances in the environment are influencing the normal working conditions of the device	
Blowing Time	2.5seconds (breath flow: $\geq 20\text{L/min}$)	
Display	128X64 black and white screen	
Dimension(LXWXT)	105mm) × 45mm × 17mm	
Batteries	7 alkaline dry batteries 1.5AAA*2	
Weight	Approx. 75g (including battery)	

7.Maintenance

Clean the surface of the equipment required, use a damp cloth or wipe a small amount of neutral detergent, avoid use of chemical organic solvents (such as alcohol, thinner, etc.) to clean the surface equipment.

Service Record

Device ID:

Purchase Date:

User:

No.	Description	Service preformed	Date
1			
2			
3			
4			
5			
6			
7			
8			

Apresys®

Apresys International Inc.

ADD: 2410 Camino Ramon, Suite 170

San Ramon, CA 94583, USA

TEL: 925.272.7788